

BreathSensor® Single Use Airflow Thermistor Now Compatible with Most PSG Systems!



The patented BreathSensor® single use Airflow Thermistor, previously compatible with only systems with a 1.5 mm AC input, such as the Sandman® SD32+™ PSG system, is now compatible with most PSG systems!

Now available on shopEmbla.com – visit now to place your order!

Features

- Available in four different sizes: Adult, Child, Infant, and Preemie
- Disposable
- Optimal for infection control
- Extremely comfortable with no nasal insert

Order Information

1. Select Your BreathSensor Size

Description	Item Number
Disposable BreathSensor - Adult	4140118
Disposable BreathSensor - Child	4140119
Disposable BreathSensor - Preemie	4140120
Disposable BreathSensor - Infant	4140121

2. Select Your BreathSensor Cable

Hardware	Item Number
S4500, N7000, Embletta®	1420446
Sandman SD32+	4140069
Sandman SD20 amplifier	4140064

For more information please contact: support@embla.com

Embla Global Headquarters • 9351 Grant St. • Suite 600 • Thornton, CO • 80229

Tel: 1 303.962.1800 • Toll-Free: 888.NO APNEA (888.662.7632) • Fax: 1.303.962.1810

Embla Amsterdam • EU Representative • Naritaweg 12 J • 1043 BZ Amsterdam • The Netherlands

Tel: +31 20 3460 130 • f: +31 20 3460 121